

产品名称: ROR beta Rabbit Monoclonal Antibody
产品货号: AMRe84265



产品概述 (Summary)

产品名称 (Production Name)	ROR beta Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ROR beta
别名 (Alternative Names)	NR1F2; ROR BETA; Rorb; RZR BETA; RZRB;; Nuclear receptor ROR beta
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q92753.A synthesized peptide derived from human Nuclear receptor ROR beta

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	53 kDa

研究背景 (Background)


产品名称: ROR beta Rabbit Monoclonal Antibody
产品货号: AMRe84265




Nuclear receptor that binds DNA as a monomer to ROR response elements (RORE) containing a single core motif half-site 5'-AGGTCA-3' preceded by a short A-T-rich sequence. Considered to have intrinsic transcriptional activity, have some natural ligands such as all-trans retinoic acid (ATRA) and other retinoids which act as inverse agonists repressing the transcriptional activity.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of ROR beta expression in HepG2 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .