

产品名称: RNaseL Rabbit Monoclonal Antibody
产品货号: AMRe84271



产品概述 (Summary)

产品名称 (Production Name)	RNaseL Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RNaseL
别名 (Alternative Names)	PRCA1; Ribonuclease4; RibonucleasL; RNase L; Rnasel; RNS4;;RNase L
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q05823.A synthesized peptide derived from human RNase L

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 84 kDa ; Observed MW: 80 kDa

研究背景 (Background)

Endoribonuclease that functions in the interferon (IFN) antiviral response. In INF treated and virus infected cells, RNASEL


产品名称: RNaseL Rabbit Monoclonal Antibody
产品货号: AMRe84271



probably mediates its antiviral effects through a combination of direct cleavage of single-stranded viral RNAs, inhibition of protein synthesis through the degradation of rRNA, induction of apoptosis, and induction of other antiviral genes.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of RNase L expression in THP1 cell lysate.

注意事项 (Note)

For research use only .