

产品名称: GEFH1 Rabbit Monoclonal Antibody
产品货号: AMRe84281



产品概述 (Summary)

产品名称 (Production Name)	GEFH1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GEFH1
别名 (Alternative Names)	ARHG2; ARHGEF2; GEF; GEFH1; Lbcl1; Lfc; LFP40;; GEF H1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q92974.A synthesized peptide derived from human GEF H1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100
蛋白分子量 (Molecular Weight)	112 kDa

研究背景 (Background)

Activates Rho-GTPases by promoting the exchange of GDP for GTP. May be involved in epithelial barrier permeability, cell


产品名称: GEFH1 Rabbit Monoclonal Antibody
产品货号: AMRe84281




motility and polarization, dendritic spine morphology, antigen presentation, leukemic cell differentiation, cell cycle regulation, innate immune response, and cancer.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

产品名称: GEFH1 Rabbit Monoclonal Antibody
产品货号: AMRe84281



注意事项 (Note)

For research use only .