

产品名称: SF3B3 Rabbit Monoclonal Antibody
产品货号: AMRe84282



产品概述 (Summary)

产品名称 (Production Name)	SF3B3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SF3B3
别名 (Alternative Names)	RSE1; SAP130; SF3b130; STAF130;;SF3B3
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q15393.A synthesized peptide derived from human SF3B3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	136 kDa

研究背景 (Background)


产品名称: SF3B3 Rabbit Monoclonal Antibody
产品货号: AMRe84282




Component of the 17S U2 SnRNP complex of the spliceosome, a large ribonucleoprotein complex that removes introns from transcribed pre-mRNAs. The 17S U2 SnRNP complex 1 directly participates in early spliceosome assembly and 2 mediates recognition of the intron branch site during pre-mRNA splicing by promoting the selection of the pre-mRNA branch-site adenosine, the nucleophile for the first step of splicing.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of SF3B3 expression in HeLa cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



产品名称: SF3B3 Rabbit Monoclonal Antibody
产品货号: AMRe84282



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .