

产品名称: KAT1/HAT1 Rabbit Monoclonal Antibody
产品货号: AMRe84284



产品概述 (Summary)

产品名称 (Production Name)	KAT1/HAT1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KAT1/HAT1
别名 (Alternative Names)	hat1; KAT1;;KAT1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O14929.A synthesized peptide derived from human KAT1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 50 kDa ; Observed MW: 45 kDa

研究背景 (Background)


产品名称: KAT1/HAT1 Rabbit Monoclonal Antibody
产品货号: AMRe84284




Histone acetyltransferase that plays a role in different biological processes including cell cycle progression, glucose metabolism, histone production or DNA damage repair. Coordinates histone production and acetylation via H4 promoter binding.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of KAT1 / HAT1 expression in (1) MCF7 cell lysate; (2) NIH/3T3 cell lysate; (3) C6 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

产品名称: KAT1/HAT1 Rabbit Monoclonal Antibody
产品货号: AMRe84284



注意事项 (Note)

For research use only .