

产品名称: USP22 Rabbit Monoclonal Antibody
产品货号: AMRe84285



产品概述 (Summary)

产品名称 (Production Name)	USP22 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	USP22
别名 (Alternative Names)	Ubiquitin carboxyl terminal hydrolase 22; Usp22; USP3L; USP22
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9UPT9.A synthesized peptide derived from human USP22

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	60 kDa

研究背景 (Background)

Deubiquitinase that plays a role in several cellular processes including transcriptional regulation, cell cycle progression or


产品名称: USP22 Rabbit Monoclonal Antibody
产品货号: AMRe84285




innate immunity. As part of the transcription regulatory histone acetylation (HAT) complex SAGA, catalyzes the deubiquitination of both histones H2A and H2B, thereby acting as a transcriptional coactivator.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of USP22 expression in (1) HeLa cell lysate; (2) Mouse spleen lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

产品名称: USP22 Rabbit Monoclonal Antibody
产品货号: AMRe84285



注意事项 (Note)

For research use only .