

产品名称: NLK Rabbit Monoclonal Antibody
产品货号: AMRe84287



产品概述 (Summary)

产品名称 (Production Name)	NLK Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NLK
别名 (Alternative Names)	LAK1; Nemo like kinase; Nlk;;NLK
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9UBE8.A synthesized peptide derived from human NLK

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, FC 1:20-1:100
蛋白分子量 (Molecular Weight)	58 kDa

研究背景 (Background)

Serine/threonine-protein kinase that regulates a number of transcription factors with key roles in cell fate


产品名称: NLK Rabbit Monoclonal Antibody
产品货号: AMRe84287



determination. Positive effector of the non-canonical Wnt signaling pathway, acting downstream of WNT5A, MAP3K7/TAK1 and HIPK2. Negative regulator of the canonical Wnt/beta-catenin signaling pathway.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of NLK expression in A375 cell lysate.

注意事项 (Note)

For research use only .