

产品名称: Retinoic Acid Receptor beta Rabbit Monoclonal Antibody
产品货号: AMRe84288

产品概述 (Summary)

产品名称 (Production Name)	Retinoic Acid Receptor beta Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Retinoic Acid Receptor beta
别名 (Alternative Names)	HAP; NR1B2; RAR B; RAR beta; RAR epsilon; RARB; RRB2;;Retinoic acid receptor beta
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P10826.A synthesized peptide derived from human Retinoic acid receptor beta

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 50 kDa ; Observed MW: 55 kDa

产品名称: Retinoic Acid Receptor beta Rabbit Monoclonal Antibody
产品货号: AMRe84288




研究背景 (Background)


Chaperonin implicated in mitochondrial protein import and macromolecular assembly. Together with Hsp10, facilitates the correct folding of imported proteins. May also prevent misfolding and promote the refolding and proper assembly of unfolded polypeptides generated under stress conditions in the mitochondrial matrix.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Retinoic Acid Receptor beta expression in MCF7 cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .