

产品名称: ICAM1/CD54 Rabbit Monoclonal Antibody
产品货号: AMRe84289



产品概述 (Summary)

产品名称 (Production Name)	ICAM1/CD54 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ICAM1/CD54 BB2; CD54; Cell surface glycoprotein P3.58; Human rhinovirus receptor;
别名 (Alternative Names)	ICAM1; intercellular adhesion molecule 1 (CD54); Intercellular adhesion molecule 1; Major group rhinovirus receptor; MALA2; MyD10; P3.58; Surface antigen of activated B cells; ICAM 1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P13597(mouse)(human:P05362).A synthesized peptide derived from mouse ICAM 2</p>

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, IHC 1:100-1:200
----------------------	-----------------------------------

产品名称: ICAM1/CD54 Rabbit Monoclonal Antibody
产品货号: AMRe84289




蛋白分子量 (Molecular Weight) Calculated MW: 58 kDa ; Observed MW: 80-110 kDa

研究背景 (Background)


ICAM proteins are ligands for the leukocyte adhesion protein LFA-1 (integrin alpha-L/beta-2). During leukocyte trans-endothelial migration, ICAM1 engagement promotes the assembly of endothelial apical cups through ARHGEF26/SGEF and RHOG activation.

研究领域 (Research Area)

图片 (Image Data)



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .