

产品名称: PLTP Rabbit Monoclonal Antibody
产品货号: AMRe84298



产品概述 (Summary)

产品名称 (Production Name)	PLTP Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PLTP
别名 (Alternative Names)	BPIFE; HDLCQ9; OD107; Pltp;;PLTP
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P55058.A synthesized peptide derived from human PLTP

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 55 kDa ; Observed MW: 80 kDa

研究背景 (Background)

Converts HDL into larger and smaller particles. May play a key role in extracellular phospholipid transport and modulation


产品名称: PLTP Rabbit Monoclonal Antibody
产品货号: AMRe84298




of hdl particles.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of PLTP expression in HeLa cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

产品名称: PLTP Rabbit Monoclonal Antibody
产品货号: AMRe84298



注意事项 (Note)

For research use only .