

产品名称: Synapsin II Rabbit Monoclonal Antibody
产品货号: AMRe84305



产品概述 (Summary)

产品名称 (Production Name)	Synapsin II Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Synapsin II
别名 (Alternative Names)	SYN 2; SYN II; SYN IIa; SYN IIb; SYN2; Synapsin2; SynapsinII; SYNII; SYNIIa; SYNIIb;;Synapsin 2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q92777.A synthesized peptide derived from human Synapsin 2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 63 kDa ; Observed MW: 75 kDa

研究背景 (Background)


产品名称: Synapsin II Rabbit Monoclonal Antibody
产品货号: AMRe84305



Neuronal phosphoprotein that coats synaptic vesicles, binds to the cytoskeleton, and is believed to function in the regulation of neurotransmitter release. May play a role in noradrenaline secretion by sympathetic neurons.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Synapsin II expression in Human fetal brain lysate.

注意事项 (Note)

For research use only .