

产品名称: LHX3 Rabbit Monoclonal Antibody
产品货号: AMRe84309



产品概述 (Summary)

产品名称 (Production Name)	LHX3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LHX3
别名 (Alternative Names)	CPHD3; LHX3; Lim3; M2 LHX3; mLim-3; mLIM3; P LIM;;LHX3
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9UBR4.A synthesized peptide derived from human LHX3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,FC 1:20-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	43 kDa

研究背景 (Background)

Transcription factor. Recognizes and binds to the consensus sequence motif 5'-AATTAATTA-3' in the regulatory elements


产品名称: LHX3 Rabbit Monoclonal Antibody
产品货号: AMRe84309



of target genes, such as glycoprotein hormones alpha chain CGA and visual system homeobox CHX10, positively modulating transcription; transcription can be co-activated by LDB2. Synergistically enhances transcription from the prolactin promoter in cooperation with POU1F1/Pit-1.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of LHX3 expression in Jurkat cell lysate.

注意事项 (Note)

For research use only .