

产品名称: TIRAP Rabbit Monoclonal Antibody
产品货号: AMRe84311



产品概述 (Summary)

产品名称 (Production Name)	TIRAP Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TIRAP
别名 (Alternative Names)	MAL; Wyatt;;TIRAP
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P58753.A synthesized peptide derived from human TIRAP

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 24 kDa ; Observed MW: 28 kDa

研究背景 (Background)

Adapter involved in TLR2, TLR4 and RAGE signaling pathways in the innate immune response. Acts via IRAK2 and TRAF-6,


产品名称: TIRAP Rabbit Monoclonal Antibody
产品货号: AMRe84311




leading to the activation of NF-kappa-B, MAPK1, MAPK3 and JNK, and resulting in cytokine secretion and the inflammatory response. Positively regulates the production of TNF-alpha (TNF) and interleukin-6 (IL6).

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of TIRAP expression in HUVEC cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

产品名称: TIRAP Rabbit Monoclonal Antibody
产品货号: AMRe84311



注意事项 (Note)

For research use only .