

产品名称: Lipin1 Rabbit Monoclonal Antibody
产品货号: AMRe84315



产品概述 (Summary)

产品名称 (Production Name)	Lipin1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Lipin1
别名 (Alternative Names)	Lipin-1; Lpin1; PAP1;;LPIN1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q14693.A synthesized peptide derived from human LPIN1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 99 kDa ; Observed MW: 130 kDa

研究背景 (Background)


产品名称: Lipin1 Rabbit Monoclonal Antibody
产品货号: AMRe84315



Acts as a magnesium-dependent phosphatidate phosphatase enzyme which catalyzes the conversion of phosphatidic acid to diacylglycerol during triglyceride, phosphatidylcholine and phosphatidylethanolamine biosynthesis and therefore controls the metabolism of fatty acids at different levels

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Lipin 1 expression in HepG2 cell lysate.

注意事项 (Note)

For research use only .