

产品名称: PICK1 Rabbit Monoclonal Antibody
产品货号: AMRe84322



产品概述 (Summary)

产品名称 (Production Name)	PICK1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PICK1
别名 (Alternative Names)	PICK 1; Pick1; PRKCABP;;PICK1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9NRD5.A synthesized peptide derived from human PICK1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,FC 1:20-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 47 kDa ; Observed MW: 50 kDa

研究背景 (Background)

Probable adapter protein that bind to and organize the subcellular localization of a variety of membrane proteins


产品名称: PICK1 Rabbit Monoclonal Antibody
产品货号: AMRe84322




containing some PDZ recognition sequence. Involved in the clustering of various receptors, possibly by acting at the receptor internalization level. Plays a role in synaptic plasticity by regulating the trafficking and internalization of AMPA receptors. May be regulated upon PRKCA activation. May regulate ASIC1/ASIC3 channel.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of PICK1 expression in HeLa cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Immunoprecipitate (IP) analysis using the Antibody at 1:50 dilution. (wb at 1:2K dilution)

产品名称: PICK1 Rabbit Monoclonal Antibody
产品货号: AMRe84322



注意事项 (Note)

For research use only .