

产品名称: OS9 Rabbit Monoclonal Antibody
产品货号: AMRe84323



产品概述 (Summary)

产品名称 (Production Name)	OS9 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	OS9
别名 (Alternative Names)	ERLEC2; Os9;;OS 9
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q13438.A synthesized peptide derived from human OS 9

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 76 kDa ; Observed MW: 83,97 kDa

研究背景 (Background)


产品名称: OS9 Rabbit Monoclonal Antibody
产品货号: AMRe84323



Lectin which functions in endoplasmic reticulum (ER) quality control and ER-associated degradation (ERAD). May bind terminally misfolded non-glycosylated proteins as well as improperly folded glycoproteins, retain them in the ER, and possibly transfer them to the ubiquitination machinery and promote their degradation. Possible targets include TRPV4.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of OS9 expression in HeLa cell lysate.

注意事项 (Note)

For research use only .