

产品名称: STUB1 Rabbit Monoclonal Antibody
产品货号: AMRe84331



产品概述 (Summary)

产品名称 (Production Name)	STUB1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, FC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	STUB1
别名 (Alternative Names)	STUB1; CHIP; HSPABP2; NY-CO-7; SDCCAG7; UBOX1;;STUB1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9UNE7.A synthesized peptide derived from human STUB1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	35 kDa

研究背景 (Background)

E3 ubiquitin-protein ligase which targets misfolded chaperone substrates towards proteasomal degradation. Plays a role in


产品名称: STUB1 Rabbit Monoclonal Antibody
产品货号: AMRe84331




the maintenance of mitochondrial morphology and promotes mitophagic removal of dysfunctional mitochondria; thereby acts as a protector against apoptosis in response to cellular stress.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of STUB1 expression in (1) HEK293 cell lysate; (2) Mouse brain lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

产品名称: STUB1 Rabbit Monoclonal Antibody
产品货号: AMRe84331



注意事项 (Note)

For research use only .