

产品名称: GRIM19 Rabbit Monoclonal Antibody
产品货号: AMRe84332



产品概述 (Summary)

产品名称 (Production Name)	GRIM19 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GRIM19
别名 (Alternative Names)	B16.6; CDA016; CGI-39; GRIM19; GRIM-19;;NDUFA13
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9P0J0.A synthesized peptide derived from human NDUFA13

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	17 kDa

研究背景 (Background)

Accessory subunit of the mitochondrial membrane respiratory chain NADH dehydrogenase (Complex I), that is believed not


产品名称: GRIM19 Rabbit Monoclonal Antibody
产品货号: AMRe84332




to be involved in catalysis. Complex I functions in the transfer of electrons from NADH to the respiratory chain. The immediate electron acceptor for the enzyme is believed to be ubiquinone.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of GRIM19 expression in Ramos cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .