

产品名称: HSP10 Rabbit Monoclonal Antibody
产品货号: AMRe84333



产品概述 (Summary)

产品名称 (Production Name)	HSP10 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HSP10
别名 (Alternative Names)	Chaperonin 10; CPN10; EPF; GROES; Hsp10; Hspe1;;Hsp10
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P61604.A synthesized peptide derived from human Hsp10

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	11 kDa

研究背景 (Background)


产品名称: HSP10 Rabbit Monoclonal Antibody
产品货号: AMRe84333



Co-chaperonin implicated in mitochondrial protein import and macromolecular assembly. Together with Hsp60, facilitates the correct folding of imported proteins. May also prevent misfolding and promote the refolding and proper assembly of unfolded polypeptides generated under stress conditions in the mitochondrial matrix.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of HSP10 expression in A431 cell lysate.

注意事项 (Note)

For research use only .