

产品名称: VEGFB Rabbit Monoclonal Antibody
产品货号: AMRe84337



产品概述 (Summary)

产品名称 (Production Name)	VEGFB Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	VEGFB
别名 (Alternative Names)	VEGF B; Vegfb; VEGFL; VRF;;VEGF B
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P49765.A synthesized peptide derived from human VEGF B

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 22 kDa ; Observed MW: 16-62 kDa

研究背景 (Background)

Growth factor for endothelial cells. VEGF-B167 binds heparin and neuropilin-1 whereas the binding to neuropilin-1 of


产品名称: VEGFB Rabbit Monoclonal Antibody
产品货号: AMRe84337



VEGF-B186 is regulated by proteolysis.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of VEGFB expression in 293T cell transfected with VEGFB recombinant protein.

注意事项 (Note)

For research use only .