

产品名称: KAISO Rabbit Monoclonal Antibody
产品货号: AMRe84341



产品概述 (Summary)

产品名称 (Production Name)	KAISO Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KAISO
别名 (Alternative Names)	Kaiso; Zbtb33; ZNF Kaiso; ZNF348;;KAISO
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q86T24.A synthesized peptide derived from human KAISO

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 74 kDa ; Observed MW: 85 kDa

研究背景 (Background)

Transcriptional regulator with bimodal DNA-binding specificity. Binds to methylated CpG dinucleotides in the consensus


产品名称: KAISO Rabbit Monoclonal Antibody
产品货号: AMRe84341



sequence 5'-CGCG-3' and also binds to the non-methylated consensus sequence 5'-CTGCNA-3' also known as the consensus kaiso binding site (KBS). Recruits the N-CoR repressor complex to promote histone deacetylation and the formation of repressive chromatin structures in target gene promoters.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of KAISO expression in MCF7 cell lysate.

注意事项 (Note)

For research use only .