

产品名称: MYF5 Rabbit Monoclonal Antibody
产品货号: AMRe84349



产品概述 (Summary)

产品名称 (Production Name)	MYF5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MYF5
别名 (Alternative Names)	Myf5; Class C basic helix loop helix protein 2;MYF5
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P13349.A synthesized peptide derived from human MYF5

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 28 kDa ; Observed MW: 37 kDa

研究背景 (Background)

Transcriptional activator that promotes transcription of muscle-specific target genes and plays a role in muscle


产品名称: MYF5 Rabbit Monoclonal Antibody
产品货号: AMRe84349



differentiation. Together with MYOG and MYOD1, co-occupies muscle-specific gene promoter core region during myogenesis. Induces fibroblasts to differentiate into myoblasts. Probable sequence specific DNA-binding protein.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of MYF5 expression in Human skeletal muscle lysate.

注意事项 (Note)

For research use only .