

产品名称: EHD1 Rabbit Monoclonal Antibody
产品货号: AMRe84351



产品概述 (Summary)

产品名称 (Production Name)	EHD1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EHD1
别名 (Alternative Names)	EHD1; H PAST; hPAST1; PAST; PAST homolog 1; PAST1; Testilin;;EHD1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9H4M9.A synthesized peptide derived from human EHD1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	61 kDa

研究背景 (Background)


产品名称: EHD1 Rabbit Monoclonal Antibody
产品货号: AMRe84351



ATP- and membrane-binding protein that controls membrane reorganization/tubulation upon ATP hydrolysis. In vitro causes vesiculation of endocytic membranes. Acts in early endocytic membrane fusion and membrane trafficking of recycling endosomes. Recruited to endosomal membranes upon nerve growth factor stimulation, indirectly regulates neurite outgrowth.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of EHD1 expression in HepG2 cell lysate.

注意事项 (Note)

For research use only .