

产品名称: IL32 Rabbit Monoclonal Antibody
产品货号: AMRe84352



产品概述 (Summary)

产品名称 (Production Name)	IL32 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IL32
别名 (Alternative Names)	IL 32alpha; IL 32beta; IL 32delta; IL 32gamma; IL32; NK4; TAIF; TAIFa; TAIFb; TAIFc; TAIFd; IL 32
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P24001.A synthesized peptide derived from human IL 32

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 27 kDa ; Observed MW: 23 kDa

研究背景 (Background)


产品名称: IL32 Rabbit Monoclonal Antibody
产品货号: AMRe84352



Cytokine that may play a role in innate and adaptive immune responses. It induces various cytokines such as TNFA/TNF-alpha and IL8. It activates typical cytokine signal pathways of NF-kappa-B and p38 MAPK.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of IL32 expression in Human tonsil lysate.

注意事项 (Note)

For research use only .