

产品名称: NUDEL Rabbit Monoclonal Antibody
产品货号: AMRe84353



产品概述 (Summary)

产品名称 (Production Name)	NUDEL Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NUDEL
别名 (Alternative Names)	EOPA; MITAP1; NDEL1; NUDEL;;NUDEL
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9GZM8.A synthesized peptide derived from human NUDEL

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 38 kDa ; Observed MW: 40 kDa

研究背景 (Background)

Required for organization of the cellular microtubule array and microtubule anchoring at the centrosome. May regulate


产品名称: NUDEL Rabbit Monoclonal Antibody
产品货号: AMRe84353




microtubule organization at least in part by targeting the microtubule severing protein KATNA1 to the centrosome. Also positively regulates the activity of the minus-end directed microtubule motor protein dynein. May enhance dynein-mediated microtubule sliding by targeting dynein to the microtubule plus ends.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of NUDEL expression in HeLa cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

产品名称: NUDEL Rabbit Monoclonal Antibody
产品货号: AMRe84353



注意事项 (Note)

For research use only .