

产品名称: CNTN2/TAG1 Rabbit Monoclonal Antibody
产品货号: AMRe84354



产品概述 (Summary)

产品名称 (Production Name)	CNTN2/TAG1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, FC
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CNTN2/TAG1
别名 (Alternative Names)	Axonin1; AXT; Cntn2; Contactin2; TAG 1; TAG564; TAX1;;Contactin 2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q02246.A synthesized peptide derived from human Contactin 2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 113 kDa ; Observed MW: 135 kDa

研究背景 (Background)

In conjunction with another transmembrane protein, CNTNAP2, contributes to the organization of axonal domains at nodes


产品名称: CNTN2/TAG1 Rabbit Monoclonal Antibody
产品货号: AMRe84354



of Ranvier by maintaining voltage-gated potassium channels at the juxtaparanodal region. May be involved in cell adhesion.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of TAG1 expression in U87-MG cell lysate.

注意事项 (Note)

For research use only .