

产品名称: Ube2N/Ubc13 Rabbit Monoclonal Antibody
产品货号: AMRe84356



产品概述 (Summary)

产品名称 (Production Name)	Ube2N/Ubc13 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Ube2N/Ubc13
别名 (Alternative Names)	BLU; HEL-S-71; Ubc13; UbCH ben; UbCH13; UBCHBEN; Ube2n;;UBE2N
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P61088.A synthesized peptide derived from human UBE2N

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	17 kDa

研究背景 (Background)

The UBE2V1-UBE2N and UBE2V2-UBE2N heterodimers catalyze the synthesis of non-canonical 'Lys-63'-linked polyubiquitin


产品名称: Ube2N/Ubc13 Rabbit Monoclonal Antibody
产品货号: AMRe84356



chains. This type of polyubiquitination does not lead to protein degradation by the proteasome. Mediates transcriptional activation of target genes. Plays a role in the control of progress through the cell cycle and differentiation.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Ube2N / Ubc13 expression in Daudi cell lysate.

注意事项 (Note)

For research use only .