

产品名称: RYK Rabbit Monoclonal Antibody
产品货号: AMRe84357



产品概述 (Summary)

产品名称 (Production Name)	RYK Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RYK
别名 (Alternative Names)	ERK 3; JTK5; JTK5A; Ryk; RYK1; Vik;;RYK
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P34925.A synthesized peptide derived from human RYK

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 68 kDa ; Observed MW: 72 kDa

研究背景 (Background)


产品名称: RYK Rabbit Monoclonal Antibody
产品货号: AMRe84357



May be a coreceptor along with FZD8 of Wnt proteins, such as WNT1, WNT3, WNT3A and WNT5A. Involved in neuron differentiation, axon guidance, corpus callosum establishment and neurite outgrowth. In response to WNT3 stimulation, receptor C-terminal cleavage occurs in its transmembrane region and allows the C-terminal intracellular product to translocate from the cytoplasm to the nucleus where it plays a crucial role in neuronal development.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of RYK expression in MCF7 cell lysate.

注意事项 (Note)

For research use only .