

产品名称: PAF1 Rabbit Monoclonal Antibody
产品货号: AMRe84360



产品概述 (Summary)

产品名称 (Production Name)	PAF1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PAF1
别名 (Alternative Names)	hPAF1; PAF1; PD2;;PAF1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q8N7H5.A synthesized peptide derived from human PAF1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 60 kDa ; Observed MW: 80 kDa

研究背景 (Background)

Component of the PAF1 complex (PAF1C) which has multiple functions during transcription by RNA polymerase II and is


产品名称: PAF1 Rabbit Monoclonal Antibody
产品货号: AMRe84360



implicated in regulation of development and maintenance of embryonic stem cell pluripotency. PAF1C associates with RNA polymerase II through interaction with POLR2A CTD non-phosphorylated and 'Ser-2'- and 'Ser-5'-phosphorylated forms and is involved in transcriptional elongation, acting both independently and synergistically with TCEA1 and in cooperation with the DSIF complex and HTATSF1.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of PAF1 expression in K562 cell lysate.

注意事项 (Note)

For research use only .