

产品名称: Bub3 Rabbit Monoclonal Antibody
产品货号: AMRe84361



产品概述 (Summary)

产品名称 (Production Name)	Bub3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Bub3
别名 (Alternative Names)	BUB3; BUB 3L;;BUB3
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O43684.A synthesized peptide derived from human BUB3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 37 kDa ; Observed MW: 40 kDa

研究背景 (Background)


产品名称: Bub3 Rabbit Monoclonal Antibody
产品货号: AMRe84361



Has a dual function in spindle-assembly checkpoint signaling and in promoting the establishment of correct kinetochore-microtubule (K-MT) attachments. Promotes the formation of stable end-on bipolar attachments. Necessary for kinetochore localization of BUB1. Regulates chromosome segregation during oocyte meiosis.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Bub3 expression in Jurkat cell lysate.

注意事项 (Note)

For research use only .