

产品名称: KIF4A Rabbit Monoclonal Antibody
产品货号: AMRe84365



产品概述 (Summary)

产品名称 (Production Name)	KIF4A Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KIF4A
别名 (Alternative Names)	Chromokinesin; KIF4; KIF4 G1; KIF4A; KIF4G1;;Chromokinesin A
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O95239.A synthesized peptide derived from human Chromokinesin A

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 140 kDa ; Observed MW: 150 kDa

研究背景 (Background)

Iron-sulfur (Fe-S) cluster binding motor protein that has a role in chromosome segregation during mitosis. Translocates


产品名称: KIF4A Rabbit Monoclonal Antibody
产品货号: AMRe84365



PRC1 to the plus ends of interdigitating spindle microtubules during the metaphase to anaphase transition, an essential step for the formation of an organized central spindle midzone and midbody and for successful cytokinesis.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of KIF4A expression in HeLa cell lysate.

注意事项 (Note)

For research use only .