

产品名称: IL23 Rabbit Monoclonal Antibody
产品货号: AMRe84369



产品概述 (Summary)

产品名称 (Production Name)	IL23 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IL23
别名 (Alternative Names)	IL 23A; IL 23p19; IL12B; IL23; IL23a; IL23P19; P19; SGRF;;IL 23A
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9NPF7.A synthesized peptide derived from human IL 23A

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	21 kDa

研究背景 (Background)

Associates with IL12B to form the pro-inflammatory cytokine IL-23 that plays different roles in innate and adaptive


产品名称: IL23 Rabbit Monoclonal Antibody
产品货号: AMRe84369



immunity. Released by antigen-presenting cells such as dendritic cells or macrophages, binds to a heterodimeric receptor complex composed of IL12RB1 and IL23R to activate JAK2 and TYK2 which then phosphorylate the receptor to form a docking site leading to the phosphorylation of STAT3 and STAT4.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of IL23 expression in Raji cell lysate.

注意事项 (Note)

For research use only .