

产品名称: FLAP Rabbit Monoclonal Antibody
产品货号: AMRe84372



产品概述 (Summary)

产品名称 (Production Name)	FLAP Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FLAP
别名 (Alternative Names)	ALOX 5AP; ALOX5 AP; ALOX5AP; FLAP;;FLAP
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P20292.A synthesized peptide derived from human FLAP

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 18 kDa ; Observed MW: 16 kDa

研究背景 (Background)

Required for leukotriene biosynthesis by ALOX5 (5-lipoxygenase). Anchors ALOX5 to the membrane. Binds arachidonic


产品名称: FLAP Rabbit Monoclonal Antibody
产品货号: AMRe84372



acid, and could play an essential role in the transfer of arachidonic acid to ALOX5. Binds to MK-886, a compound that blocks the biosynthesis of leukotrienes.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of FLAP expression in THP1 cell lysate.

注意事项 (Note)

For research use only .