

产品名称: p60 CAF1 Rabbit Monoclonal Antibody
产品货号: AMRe84378



产品概述 (Summary)

产品名称 (Production Name)	p60 CAF1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	p60 CAF1
别名 (Alternative Names)	CAF 1A; CAF IP60; CAF-1; CAF1; CAF1A; CAF1P60; CHAF1B; MPHOSPH7; MPP7;;CAF 1 subunit B
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q13112.A synthesized peptide derived from human CAF 1 subunit B

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 61 kDa ; Observed MW: 70 kDa

产品名称: p60 CAF1 Rabbit Monoclonal Antibody
产品货号: AMRe84378




研究背景 (Background)


Complex that is thought to mediate chromatin assembly in DNA replication and DNA repair. Assembles histone octamers onto replicating DNA in vitro. CAF-1 performs the first step of the nucleosome assembly process, bringing newly synthesized histones H3 and H4 to replicating DNA; histones H2A/H2B can bind to this chromatin precursor subsequent to DNA replication to complete the histone octamer.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of p60 CAF1 expression in Caco2 cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .