

产品名称: NAIP Rabbit Monoclonal Antibody
产品货号: AMRe84380



产品概述 (Summary)

产品名称 (Production Name)	NAIP Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NAIP
别名 (Alternative Names)	BIRC1; Birc1a; NAIP; Naip1; NLRB1; psiNAIP;;NAIP
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q13075.A synthesized peptide derived from human NAIP

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	160 kDa

研究背景 (Background)

Prevents motor-neuron apoptosis induced by a variety of signals. Possible role in the prevention of spinal muscular atrophy


产品名称: NAIP Rabbit Monoclonal Antibody
产品货号: AMRe84380



that seems to be caused by inappropriate persistence of motor-neuron apoptosis: mutated or deleted forms of NAIP have been found in individuals with severe spinal muscular atrophy.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of NAIP expression in Human fetal kidney lysate.

注意事项 (Note)

For research use only .