

产品名称: GRIK2 Rabbit Monoclonal Antibody
产品货号: AMRe84386



产品概述 (Summary)

产品名称 (Production Name)	GRIK2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GRIK2
别名 (Alternative Names)	EAA4; GLR6; MRT6; GLUK6; GLUR6; GluK2;;GluK2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q13002.A synthesized peptide derived from human GluK2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	103 kDa

研究背景 (Background)

Ionotropic glutamate receptor. L-glutamate acts as an excitatory neurotransmitter at many synapses in the central nervous


产品名称: GRIK2 Rabbit Monoclonal Antibody
产品货号: AMRe84386



system. Binding of the excitatory neurotransmitter L-glutamate induces a conformation change, leading to the opening of the cation channel, and thereby converts the chemical signal to an electrical impulse.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of GRIK2 expression in A431 cell lysate.

注意事项 (Note)

For research use only .