

产品名称: Monoacylglycerol Lipase Rabbit Monoclonal Antibody
产品货号: AMRe84387



产品概述 (Summary)

产品名称 (Production Name)	Monoacylglycerol Lipase Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Monoacylglycerol Lipase
别名 (Alternative Names)	HUK5; Lysophospholipase homolog; Lysophospholipase like; MAGL; MGL; MGLL; Monoacylglycerol lipase;; Monoglyceride lipase
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q99685.A synthesized peptide derived from human Monoglyceride lipase

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	33 kDa

产品名称: Monoacylglycerol Lipase Rabbit Monoclonal Antibody
产品货号: AMRe84387




研究背景 (Background)

Converts monoacylglycerides to free fatty acids and glycerol. Hydrolyzes the endocannabinoid 2-arachidonoylglycerol, and thereby contributes to the regulation of endocannabinoid signaling, nociception and perception of pain.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Monoacylglycerol Lipase expression in human heart lysate.

注意事项 (Note)

For research use only .