

产品名称: ZNF217 Rabbit Monoclonal Antibody
产品货号: AMRe84388



产品概述 (Summary)

产品名称 (Production Name)	ZNF217 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ZNF217
别名 (Alternative Names)	ZABC1; ZN217; ZNF217;;ZNF217
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O75362.A synthesized peptide derived from human ZNF217

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 115 kDa ; Observed MW: 125 kDa

研究背景 (Background)

Binds to the promoters of target genes and functions as repressor. Promotes cell proliferation and antagonizes cell death.


产品名称: ZNF217 Rabbit Monoclonal Antibody
产品货号: AMRe84388



Promotes phosphorylation of AKT1 at 'Ser-473'.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of ZNF217 expression in Caco2 cell lysate.

注意事项 (Note)

For research use only .