

产品名称: c-Maf Rabbit Monoclonal Antibody
产品货号: AMRe84396



产品概述 (Summary)

产品名称 (Production Name)	c-Maf Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	c-Maf
别名 (Alternative Names)	cMaf; maf; MAF2; Proto oncogene c Maf;;c Maf
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O75444.A synthesized peptide derived from human c Maf

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 38 kDa ; Observed MW: 44 kDa

研究背景 (Background)

Acts as a transcriptional activator or repressor. Involved in embryonic lens fiber cell development. Recruits the


产品名称: c-Maf Rabbit Monoclonal Antibody
产品货号: AMRe84396




transcriptional coactivators CREBBP and/or EP300 to crystallin promoters leading to up-regulation of crystallin gene during lens fiber cell differentiation.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of c-Maf expression in (1) HUVEC cell lysate; (2) NIH/3T3 cell lysate; (3) PC12 cell lysate.




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

产品名称: c-Maf Rabbit Monoclonal Antibody
产品货号: AMRe84396



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .