

产品名称: Methionine Aminopeptidase 2 Rabbit Monoclonal Antibody
产品货号: AMRe84402



产品概述 (Summary)

产品名称 (Production Name)	Methionine Aminopeptidase 2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Methionine Aminopeptidase 2
别名 (Alternative Names)	Amp2; MAP2; Metap2; MNPEP; p67; p67eIF2;;METAP2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P50579.A synthesized peptide derived from human METAP2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 53 kDa ; Observed MW: 67 kDa

产品名称: Methionine Aminopeptidase 2 Rabbit Monoclonal Antibody
产品货号: AMRe84402




研究背景 (Background)


Cotranslationally removes the N-terminal methionine from nascent proteins. The N-terminal methionine is often cleaved when the second residue in the primary sequence is small and uncharged (Met-Ala-, Cys, Gly, Pro, Ser, Thr, or Val). The catalytic activity of human METAP2 toward Met-Val peptides is consistently two orders of magnitude higher than that of METAP1, suggesting that it is responsible for processing proteins containing N-terminal Met-Val and Met-Thr sequences in vivo.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Methionine Aminopeptidase 2 expression in HeLa cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .