

产品名称: PCK1 Rabbit Monoclonal Antibody
产品货号: AMRe84405



产品概述 (Summary)

产品名称 (Production Name)	PCK1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PCK1
别名 (Alternative Names)	GTP; PCK1; PEP carboxykinase; PEPCK1; PEPCKC;;PCK1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P35558.A synthesized peptide derived from human PCK1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	69 kDa

研究背景 (Background)

Catalyzes the conversion of oxaloacetate (OAA) to phosphoenolpyruvate (PEP), the rate-limiting step in the metabolic


产品名称: PCK1 Rabbit Monoclonal Antibody
产品货号: AMRe84405



pathway that produces glucose from lactate and other precursors derived from the citric acid cycle.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of PCK1 expression in Human fetal kidney lysate.

注意事项 (Note)

For research use only .