

产品名称: MADCAM1 Rabbit Monoclonal Antibody
产品货号: AMRe84407



产品概述 (Summary)

产品名称 (Production Name)	MADCAM1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MADCAM1
别名 (Alternative Names)	MACAM1; Madcam1;;MADCAM1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q13477.A synthesized peptide derived from human MADCAM1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 40 kDa ; Observed MW: 60 kDa

研究背景 (Background)

Cell adhesion leukocyte receptor expressed by mucosal venules, helps to direct lymphocyte traffic into mucosal tissues


产品名称: MADCAM1 Rabbit Monoclonal Antibody
产品货号: AMRe84407



including the Peyer patches and the intestinal lamina propria. It can bind both integrin alpha-4/beta-7 and L-selectin, regulating both the passage and retention of leukocytes. Isoform 2, lacking the mucin-like domain, may be specialized in supporting integrin alpha-4/beta-7-dependent adhesion strengthening, independent of L-selectin binding.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of MADCAM1 expression in HUVEC cell lysate.

注意事项 (Note)

For research use only .