

产品名称: beta Defensin 1 Rabbit Monoclonal Antibody
产品货号: AMRe84409



产品概述 (Summary)

产品名称 (Production Name)	beta Defensin 1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	beta Defensin 1
别名 (Alternative Names)	BD1; beta 1; DEFB1; DEFB101; Defensin; Defensin beta 1; HBD1;;beta Defensin 1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P60022.A synthesized peptide derived from human beta Defensin 1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	7 kDa

研究背景 (Background)


产品名称: beta Defensin 1 Rabbit Monoclonal Antibody
产品货号: AMRe84409



Has bactericidal activity. May act as a ligand for C-C chemokine receptor CCR6. Positively regulates the sperm motility and bactericidal activity in a CCR6-dependent manner. Binds to CCR6 and triggers Ca^{2+} mobilization in the sperm which is important for its motility.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of beta Defensin 1 expression in Human recombinant beta Defensin 1 lysate.

注意事项 (Note)

For research use only .