

产品名称: CLPP Rabbit Monoclonal Antibody
产品货号: AMRe84414



产品概述 (Summary)

产品名称 (Production Name)	CLPP Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CLPP
别名 (Alternative Names)	Endopeptidase Clp;;CLPP
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q16740.A synthesized peptide derived from human CLPP

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 30 kDa ; Observed MW: 26 kDa

研究背景 (Background)

Protease component of the Clp complex that cleaves peptides and various proteins in an ATP-dependent process. Has low


产品名称: CLPP Rabbit Monoclonal Antibody
产品货号: AMRe84414




peptidase activity in the absence of CLPX. The Clp complex can degrade CSN1S1, CSN2 and CSN3, as well as synthetic peptides (in vitro) and may be responsible for a fairly general and central housekeeping function rather than for the degradation of specific substrates.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of CLPP expression in A431 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .