

产品名称: FABPI Rabbit Monoclonal Antibody
产品货号: AMRe84415



产品概述 (Summary)

产品名称 (Production Name)	FABPI Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FABPI
别名 (Alternative Names)	FABP2; FABPI; IFABP;;Fatty acid binding protein 2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P12104.A synthesized peptide derived from human Fatty acid binding protein 2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	15 kDa

研究背景 (Background)


产品名称: FABPI Rabbit Monoclonal Antibody
产品货号: AMRe84415



FABPs are thought to play a role in the intracellular transport of long-chain fatty acids and their acyl-CoA esters. FABP2 is probably involved in triglyceride-rich lipoprotein synthesis. Binds saturated long-chain fatty acids with a high affinity, but binds with a lower affinity to unsaturated long-chain fatty acids.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of intestinal FABP expression in Human small intestine lysate.

注意事项 (Note)

For research use only .