

产品名称: TAF1 Rabbit Monoclonal Antibody
产品货号: AMRe84416



产品概述 (Summary)

产品名称 (Production Name)	TAF1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TAF1
别名 (Alternative Names)	BA2R; CCG1; CCGS; DYT3; N TAF1; NSCL2; p250; TAF1; TAF2A; TAFII250; XDP;;TAF1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P21675.A synthesized peptide derived from human TAF1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 213 kDa ; Observed MW: 300 kDa

研究背景 (Background)


产品名称: TAF1 Rabbit Monoclonal Antibody
产品货号: AMRe84416



Largest component and core scaffold of the TFIID basal transcription factor complex. Contains novel N- and C-terminal Ser/Thr kinase domains which can autophosphorylate or transphosphorylate other transcription factors.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of TAF1 expression in SH-SY5Y cell lysate.

注意事项 (Note)

For research use only .