

产品名称: SNX1 Rabbit Monoclonal Antibody
产品货号: AMRe84417



产品概述 (Summary)

产品名称 (Production Name)	SNX1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SNX1
别名 (Alternative Names)	Snx1; SNX1A; Sorting nexin 1; Sorting nexin 1A; Vps5;;Sorting nexin 1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q13596.A synthesized peptide derived from human Sorting nexin 1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 59 kDa ; Observed MW: 74 kDa

研究背景 (Background)


产品名称: SNX1 Rabbit Monoclonal Antibody
产品货号: AMRe84417



May be involved in several stages of intracellular trafficking. Plays a role in targeting ligand-activated EGFR to the lysosomes for degradation after endocytosis from the cell surface and release from the Golgi.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of SNX1 expression in 293T cell lysate.

注意事项 (Note)

For research use only .