

产品名称: Cytochrome P450 2C9 Rabbit Monoclonal Antibody
产品货号: AMRe84424



产品概述 (Summary)

产品名称 (Production Name)	Cytochrome P450 2C9 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Cytochrome P450 2C9
别名 (Alternative Names)	CPC9; CYP2C; CYP2C10; CYP2C9; CYP11C9; P450 PB1; P450MP;;Cytochrome P450 2C9
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P11712.A synthesized peptide derived from human Cytochrome P450 2C9

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 56 kDa ; Observed MW: 50 kDa

产品名称: Cytochrome P450 2C9 Rabbit Monoclonal Antibody
产品货号: AMRe84424




研究背景 (Background)

A cytochrome P450 monooxygenase involved in the metabolism of various endogenous substrates, including fatty acids and steroids.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Cytochrome P450 2C9 expression in HepG2 cell lysate.

注意事项 (Note)

For research use only .