

产品名称: ACVR2A Rabbit Monoclonal Antibody
产品货号: AMRe84430



产品概述 (Summary)

产品名称 (Production Name)	ACVR2A Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ACVR2A
别名 (Alternative Names)	ACTR2; ACTRII; ACTRIIA; Acvr2; ACVR2A;;ACVR2A
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P27037.A synthesized peptide derived from human ACVR2A

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 58 kDa ; Observed MW: 85 kDa

研究背景 (Background)

On ligand binding, forms a receptor complex consisting of two type II and two type I transmembrane serine/threonine


产品名称: ACVR2A Rabbit Monoclonal Antibody
产品货号: AMRe84430




kinases. Type II receptors phosphorylate and activate type I receptors which autophosphorylate, then bind and activate SMAD transcriptional regulators. Receptor for activin A, activin B and inhibin A.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

产品名称: ACVR2A Rabbit Monoclonal Antibody
产品货号: AMRe84430



注意事项 (Note)

For research use only .