

产品名称: MYO7A Rabbit Monoclonal Antibody
产品货号: AMRe84433



产品概述 (Summary)

产品名称 (Production Name)	MYO7A Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MYO7A
别名 (Alternative Names)	DFNA11; DFNB2; Myo7a; MYU7A; NSRD2; Ush1b;; Myosin VIIa
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q13402.A synthesized peptide derived from human Myosin VIIa

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100
蛋白分子量 (Molecular Weight)	254 kDa

研究背景 (Background)

Myosins are actin-based motor molecules with ATPase activity. Unconventional myosins serve in intracellular movements.


产品名称: MYO7A Rabbit Monoclonal Antibody
产品货号: AMRe84433



Their highly divergent tails are presumed to bind to membranous compartments, which would be moved relative to actin filaments.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Myosin VIIa expression in Y79 cell lysate.

注意事项 (Note)

For research use only .